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
Chart projector

Mod. CP 2047, CP 2147, CP 2033,

CP 3000, CP 3137, CP 1042, CP 1142

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CP 2047

Two disks. Projection enlargement x30 at 5 mt. 45 tests, masks of horizontal, vertical lines and single letters, filters that can be superimposed on all the tests: red/green, polarising filters. Optotypes: Sloan's letters (according to the UNI EN ISO 8596-8597 specifications), Landolt's rings, E directionals, numbers, drawings for children. LogMAR progression (according to the UNI EN ISO 8596-8597 specifications), decimal notation.

Mod. CP 2047

The chart displays 45 tests arranged in a grid. The top two rows contain Sloan letters and Landolt rings. The middle section includes directional 'E's, numbers, and child drawings. The bottom section features various filters and alignment tests. Each test is accompanied by its logMAR value.

E .05	C .08	D V H .1 C N R .125	N C K Z O .16 R H S D K .2 C Z R H S .25	O N H R C .32 D K S N V .4 Z S O K N .5	H Z O V C .63 C R H D V .8 Z K C S N 1.0	V H C N O 1.0 S V Z D K 1.25 D S R K N 1.6	2 5 8 .1 7 3 4 .15	9 0 2 6 8 .2 3 7 5 0 4 .3 2 8 9 3 5 .4
7 0 4 2 8 .5 3 6 5 7 9 .6 0 2 3 8 4 .7	4 0 5 7 3 .8 6 2 9 8 0 .9 3 7 4 5 2 1.0	E .05	W .08	E W M .1 E W M .125	W E E M E .16 M W E E W .2 E M W M E .25	W M W M E .32 W M E E W .4 E M E E W .5	W E E M W .63 M E E W M .8 E M W M E 1.0	C .08
C O O .1 O C O .125	O O O O O .16 O O O O O .2 O O O O O .25	O O O O O .32 O O O O O .4 O O O O O .5	O O O O O .63 O O O O O .8 O O O O O 1.0	O O O O O 1.25 O O O O O 1.6 O O O O O 2.0	Child drawings .1	Child drawings .2	Child drawings .4	Child drawings .7
Duochrome balance	Coincidence vertical	Vertical line	Horizontal line	Red/green	Phoria	Phoria with fixation	Worth	Minute stereo
Astigmatism clock dial	Coincidence horizontal	Cross cylinder	Fixation disparity	Schober	Binocular balance	Cross grid	Schematic eye	Fixation

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CP 2147

Two disks. Projection enlargement x30 at 5 mt. 45 tests, masks of horizontal, vertical lines and single letters, filters that can be superimposed on all the tests: red/green, polarising filters. Optotypes: Sloan's letters (according to the UNI EN ISO 8596-8597 specifications), Landolt's rings, E directionals, numbers, drawings for children. Linear progression of the optotypes.

Mod. CP 2147

The chart displays 45 tests arranged in a grid. The top two rows contain Sloan letters and Landolt rings. The middle section includes directional 'E's, numbers, and child drawings. The bottom section features various filters and alignment tests. Each test is accompanied by its logMAR value.


E .05	C .08	D H .1 N .1	N K O .15 R H S D K .2 C Z R H S .3	O N H R C .4 D K S N V .5 Z S O K N .6	H Z O V C .7 C R H D V .8 Z K C S N .9	V H C N O 1.0 S V Z D K 1.2 D S R K N 1.5	2 5 8 .1 7 3 4 .15	9 0 2 6 8 .2 3 7 5 0 4 .3 2 8 9 3 5 .4
7 0 4 2 8 .5 3 6 5 7 9 .6 0 2 3 8 4 .7	4 0 5 7 3 .8 6 2 9 8 0 .9 3 7 4 5 2 1.0	E .05	W .1	E W M .15 E W M .2	W E E M E .3 M W E E W .4 E M W M E .5	W M W M E .6 W M E E W .7 E M E E W .8	W E E M W .9 M E E W M 1.0 E M W M E 1.25	C .08
C O O .1 O C O .1	O O O O O .15 O O O O O .2 O O O O O .3	O O O O O .4 O O O O O .5 O O O O O .6	O O O O O .7 O O O O O .8 O O O O O .9	O O O O O 1.0 O O O O O 1.5 O O O O O 2.0	Child drawings .1	Child drawings .2	Child drawings .4	Child drawings .7
Duochrome balance	Coincidence vertical	Vertical line	Horizontal line	Red/green	Phoria	Phoria with fixation	Worth	Minute stereo
Astigmatism clock dial	Coincidence horizontal	Cross cylinder	Fixation disparity	Schober	Binocular balance	Cross grid	Schematic eye	Fixation

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CP 2033 ZOOM (USA model)

Two disks. Zoom projection enlargement. 33 tests, masks of horizontal, vertical lines and single letters, filters that can be superimposed on all the tests: red/green, polarising filters. Optotypes: Sloan's letters (according to the ANSI specifications), E directional, numbers, drawings for children. Progression according to Snellen notation (20 ft).

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
CP 1042

One disk. Projection enlargement x30 at 5 mt. 42 tests, no superimposable filter or masking. Optotypes: Sloan's letters (according to the UNI EN ISO 8596 - 8597 specifications), Landolt's rings, E directional, numbers, drawings for children. LogMAR progression (according to the UNI EN ISO 8596 - 8597 specifications), decimal notation.

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CP 1142

One disk. Projection enlargement x30 at 5 mt. 42 tests, no superimposable filter or masking. Optotypes: Sloan's letters (according to the UNI EN ISO 8596 - 8597 specifications), Landolt's rings, E directionals, numbers, drawings for children. Linear progression of the optotypes.

Mod. CP 1142

0.05 E 0.08 C 0.1 D H 0.15 N K O 0.2 RHSDK 0.3 CZRHS 0.4 ONHRC 0.5 DKSNV 0.6 ZSOKN 0.7 HZOV 0.8 CRHDV 0.9 ZKCSN 1.0 VHCNO 1.2 SVZDK 1.5 DSRKN 1.1 2 5 8 1.2 9 0 2 6 8 1.3 3 7 5 0 4 1.4 7 3 4 1.5 2 8 9 3 5 4

0.5 7 0 4 2 8 0.6 3 6 5 7 9 0.7 0 2 3 8 4 0.8 4 0 5 7 3 0.9 6 2 9 8 0 1.0 3 7 4 5 2

0.05 E 0.1 W 0.15 EWM 0.2 EEW 0.3 WEEME 0.4 MWEEM 0.5 EMMEM 0.6 EMMEM 0.7 WEMEM 0.8 MEMEM 0.9 EMMEM 1.0 WEMEM 1.25 C 0.8 C 0.1 C 0.1 C 0.15 C C C C 0.2 C C C C 0.3 C C C C 0.4 C C C C 0.5 C C C C 0.6 C C C C 0.7 C C C C 0.8 C C C C 0.9 C C C C 1.0 C C C C 1.5 C C C C 2.0 C C C C

0.1 C C 0.1 C C 0.15 C C C C 0.2 C C C C 0.3 C C C C 0.4 C C C C 0.5 C C C C 0.6 C C C C 0.7 C C C C 0.8 C C C C 0.9 C C C C 1.0 C C C C 1.5 C C C C 2.0 C C C C

0.1 6 0 8 8 0 5 0.5 8 0 5 5 0 3 0.15 O N H R C 0.2 R H S D K 0.3 C Z R H S 0.4 6 0 8 8 0 5 0.5 8 7 5 2 2 5 7 0.1 Phoria with fixation 0.2 Worth 0.3 Minute stereo

Duochrome balance Coincidence vertical Vertical/horizontal line Red/green Phoria with fixation Worth Minute stereo

0.1 Astigmatism clock dial 0.1 Coincidence horizontal 0.15 Cross cylinder 0.2 Fixation disparity 0.3 Schober 0.4 Binocular balance 0.5 Cross grid 0.6 Schematic eye

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CP 3000

Two disks. Projection enlargement x23,89 at 5 mt. 37 tests, masks of horizontal lines and single letters, filters that can be superimposed on all the tests: red/green, polarising filters. Optotypes: Sloan's letters (according to the UNI EN ISO 8596-8597 specifications), Landolt's rings, E directionals, numbers, drawings for children. LogMAR progression (according to the UNI EN ISO 8596-8597 specifications), decimal notation.

Mod. CP 3000

400 E 0.05 200 FZ 0.1 25 LPH 0.16 80 VHRN 0.25 40 LNEO 0.5 20 ERPD 1.0 100 0.2 40 0.5 25 0.8 160 D N 0.125 100 FVOR 0.2 50 ECDH 0.4 25 0.25 0.3 0.6 0.9 1.0 1.25 1.6 1.25 1.6 1.0 0.25 0.4 0.5 0.63 0.7 1.0 1.25 1.6 2.0

200 2 6 0.1 70 9 7 8 3 0.3 30 4 6 0 9 0.6 22 6 5 2 4 0.9 160 C C 0.125 100 C C C C 0.2 50 C C C C 0.4 28 C C C C 0.7 16 C C C C 1.25 100 8 2 7 0.2 40 0 2 4 6 0.4 28 3 7 8 2 0.7 20 2 4 6 8 1.0 16 5 7 9 3 1.25 125 C C C 0.16 60 C C C C 0.25 40 C C C C 0.5 25 C C C C 0.8 12.5 C C C C 0.63 30 C C C C 1.0 20 C C C C 1.6

125 EWM 0.16 80 WEMEM 0.25 40 MEMEM 0.5 20 MEMEM 1.0 0.8 8 0 8 0 5 0.5 7 5 2 2 5 7 0.1 Astigmatism clock dial 0.2 Binocular balance 0.3 Worth 0.4 Minute stereo

100 EEW 0.2 50 WEMEM 0.4 25 MEMEM 0.8 12.5 EWMEM 1.6 0.1 Cross grid 0.2 Schober 0.3 Cross cylinder 0.4 Duochrome balance 0.5 Phoria 0.6 Phoria with fixation 0.7 Coincidence Horizontal 0.8 Coincidence vertical 0.9 Fixation 1.0 Vertical line

24	INTRODUCTION
	STANDARD EQUIPMENT
25	DESCRIPTION AND USE
	AMBIENT CONDITIONS FOR SHIPPING, STORAGE, AND OPERATION
26	END-OF-LIFE DISPOSAL INFORMATION
27	LEGENDA
28	INSTALLATION AND COMMISSIONING
	Installing the Projector
29	Installing Multiple Projectors in a Single Room
	Positioning the Projector and Focusing the Image on the Screen
30	SAFETY PRECAUTIONS, WARNINGS, AND GENERAL SAFETY MEASURES
31	INSTRUCTIONS FOR USE
	Programming Test Sequences
	CP 3000 and CP 3137 Projectors
	CP 2047, CP 2147, CP 1042, CP 1142, and CP 2033 Projectors
32	Projection Distance
	Projection with Mirrors
33	OPERATIONAL FEATURES
	Variable-Contrast Chart Projector (CP 3000 and CP 3137)
	Dual-Disk Chart Projector, 30x enlargement (CP 2047 and CP 2147)
	Single-Disk Chart Projector, 30x enlargement (CP 1042 and CP 1142)
	Zoom Chart Projector (CP 2033)
	TECHNICAL SPECIFICATIONS
	Classification According to EN60601-1 Standard
34	Block Diagram of Operation
35	ROUTINE MAINTENANCE
	Replacing the Lamp
	Replacing the Remote Control Batteries
	Replacing the Line Fuses
	Protecting the Projector from Dust, Cleaning
	Other Precautionary Measures
36	LIABILITY
	WARRANTY AND TECHNICAL ASSISTANCE
37	EXPLANATION OF GRAPHICS SYMBOLS
38	STATEMENTS OF COMPLIANCE
39	REFERENCE STANDARDS
40	CERTIFICATIONS

READ THIS MANUAL
CAREFULLY

We congratulate you on your purchase of the CP Chart Projector by CSO


All our products are manufactured with maximum attention to safety. The chart projector is a high-performance instrument. In order to use it efficiently and in complete safety we recommend reading this manual carefully before installation and use and heeding all the safety warnings contained herein and on the instrument labels. Even operators who are familiar with this type of instrument are urged to reread the instructions contained in this manual. Always keep the manual in proximity to the instrument for consultation during use.


This document constitutes the complete manual.

STANDARD EQUIPMENT


On receipt of the instrument, remove it from its packing and check that all the components listed below are included:

- a) Projector, complete with accessories for requested installation.
- b) One infrared remote control (13).
- c) One pair red/green viewing glasses (14).
- d) One pair polarized viewing glasses (15).
- e) One power cord.
- f) This manual.
- g) One wall-mount projection screen (17).
- h) One dust cover (18).
- i) One T wrench for adjusting the projection angle (19).
- j) One socket wrench for opening the casing (20).
- k) One spare halogen lamp, 12V 20W G4 (10), two spare halogen lamps 6V 25W G4 (21) for the contrast projector model.
- l) Two spare fuses (22).
- m) 4 batteries for the remote control, 1.5V AAA (21).

Data plate on models CP 2047,
CP 2147, 1042, 1142, and CP 2033

Data plate on models CP 3000 and
CP 3137

DESCRIPTION AND USE

The instruments in the CHART PROJECTOR family are designed for projection of optotypes on a wall-mounted projection screen. Operation of the instrument is analogous to that of traditional chart projectors. The instruments can project high-quality optotypes, created using sophisticated DIADUR technology, for use by the ophthalmologist or optician, each within the limits of his competence, in making measurements and clinical evaluations of visual acuity and point up limitations and defects in vision.

This type of instrument permits projecting up to 47 different tests (according to model) including stereoscopic testing, tests for heterophoria and astigmatism, red-green tests, and polarized tests, in addition to letters, numbers, directional Es, figures for children, and Landolt rings.

The instruments are designed to be used only in spaces for medical use (in particular, spaces suitable for use of Type "B" devices with applied parts); for example, in ophthalmology clinics or in optician's or optometry.

laboratories and/or centers. The electrical system powering the instrument must comply with the relevant laws and regulations in force in the country of use.

AMBIENT CONDITIONS FOR SHIPPING, STORAGE, AND OPERATION

As long as the instrument remains in its original packing it may be exposed to the environmental conditions listed below, for a maximum of 15 weeks during shipping and warehousing, without suffering damage:

Temperature between -10 °C and +60 °C;
atmospheric pressure between 500 hPa and 1060 hPa;
relative humidity between 10% and 90%.

Ambient conditions for operation are:

Temperature between +15 °C and +30 °C;
atmospheric pressure between 700 hPa and 1060 hPa;
relative humidity between 30% and 75%.

The symbol (barred waste can) shown in the figure and found on the exterior of the instrument indicates that the electrical and electronic parts of the end-of-life instrument must be **separated and disposed of as special waste.**

END-OF-LIFE DISPOSAL INFORMATION

in accordance with Art. 13 of Leg. Decree No. 151 of 25 July 2005 implementing Directives 2002/95/EC, 2002/96/EC, and 2003/108/EC concerning reduction of use of dangerous substances in electrical and electronic equipment and disposal of electrical and electronic waste.

The instrument you have purchased is made using particular materials and substances. It may also contain substances having potentially dangerous effects on the environment and human health if released into the environment by improper disposal.

To avoid releasing dangerous substances into the environment and in order to promote conservation of natural resources, and should the user decide to dispose of an end-of-life instrument, the manufacturer will facilitate re-use and recovery and recycling of the materials it contains.

Government agencies have adopted measures obliging users, distributors, and manufacturers to contribute to collection of waste electrical and electronic equipment (WEEE) and prescribes that such equipment find re-use or be recovered or recycled.

When disposing of the instrument, remember that disposal is regulated by precise European and national laws and regulations that prescribe the following:

- **Do not dispose of as ordinary municipal waste.** For separate disposal contact a company specialized in disposal of waste electrical and electronic equipment or your local waste disposal agency for information.

- If a new instrument is purchased from the same manufacturer to replace an end-of-life instrument put on the market before 13 August 2005, of an equivalent type and performing the same functions as the new device, the distributor or manufacturer is required by law to take back the end-of-life device.




- If the user intends to dispose of a used device put on the market after 13 August 2005, the distributor or manufacturer is required by law to take back the device.

- The manufacturer is responsible for transporting, treating, and recovering and/or recycling any used equipment collected and for all relevant expenses.

- **Never forget that dangerous substances present in waste electrical and electronic equipment and/or improper use of same or parts of same can have potentially adverse effects on the environment and human health.** The instrument described in this manual contains metal and plastic mechanical parts, electrical components, and electronic circuit cards. The manufacturer is at user's complete disposition for any information requested regarding the dangerous substances contained in the instrument and recovery and recycling procedures and/or the possibility of re-using the end-of-life instrument.

Current legislation provides severe sanctions in the case of failure to respect disposal laws and regulations in force.

- 1) Casing
- 2) Infrared receiver
- 3) Projector lens
- 4) Focusing ring
- 5) Casing screws
- 6) Tilt adjustment lock screw
- 7) Lighted power switch
- 8) Power socket with fuse/voltage selector unit
- 9) Data plate
- 10) Halogen lamp, 12V 20W G4, code no. 3008011220C
- 11) Lamp cover
- 12) Halogen lamp, 6V 25W G4, code no. 3008010625A
- 13) Remote control
- 14) Red-green viewing glasses, code no. 100301050
- 15) Polarized viewing glasses, code no. 100301049
- 16) Contrast projector display
- 17) Projection screen, code no. 100301048
- 18) Dust cover, code no. 4013110
- 19) T wrench for tilt adjustment, code no. 6015T6
- 20) Socket wrench, 3
- 21) Batteries, 1.5V AAA
- 22) Spare fuses
- 23) RS 232 serial interface socket (not all models)
- 24) Focus adjustment knob on Zoom projector (CP 2033)
- 25) Projection enlargement knob on Zoom projector (CP 2033)


Attention

Read this manual, and this chapter in particular, carefully before switching on and using the instrument.

The instrument is supplied packed in such a manner as to best withstand standard shipping and warehousing conditions. Should you notice defects attributable to shipping when unpacking the instrument, contact your installation service or the manufacturer.

Installing the Projector


The instrument is supplied completely assembled and ready for the type of installation requested:

- **Wall mounting:** secure the plate to the wall with screws. Install the projector support arm (similar to that used on refraction units) to the plate.
- **Mounting on refraction unit:** insert the mounting pin and lock in place using the wrench supplied with the instrument; insert the pin in the support arm on the refraction unit.
- **Table-top installation:** insert the sphere and lock in place using the wrench supplied with the instrument.

Whatever the type of installation configuration ordered, the vertical and/or horizontal tilt of the projector will always require adjustment. Use the T wrench (19) supplied with the projector to loosen the screws (6). Bring the projector into perfect alignment with the projection screen (17) and retighten the screws (6).

Check that the voltage selector on the power socket (8) is set to the voltage at which the instrument will be used (120V AC or 230 V AC). If necessary, change the voltage setting as explained below:

- Extract the fuse-holder box.
- Rotate until the desired voltage value appears in the window.
- Check that the fuse rating is compatible with working voltage (see TECHNICAL SPECIFICATIONS) and if necessary replace the fuses.
- Replace the box.


Voltage selector in the power socket.

Check that the projector is switched off; plug the power cord into the line socket.

- Check that your line voltage is the same as that reported on the instrument data plate. Should the values differ, contact your customer service or the manufacturer. Your entire electrical system must comply with the CEI 64-4 standard or the more recent CEI 64-8 sect. 710 (Electrical Systems for Medical Use). In case of doubt, contact your electrical system installation and maintenance service.
- Never use multiple plugs, adapters, or extension cords to connect the instrument plug to the line socket.
- When unplugging the instrument from line power supply, even in an emergency, grasp the plug only; never pull the cord to disconnect the plug.

Installing Multiple Projectors in a Single Room

All our projectors are programmed for remote control transmission using the code 00; in other words, the default transmission code is the same for all remotes. In the case more than one projector are used in a single room, the remote control transmission codes must be modified to avoid interference among them. Open the battery compartment on the rear of the remote (refer to figure on left) and enter a new combination by resetting one DIP switch at a time to ON (so obtaining four additional combinations beside the standard code). When finished, turn the projector off (main switch, 7). Bring the remote to a distance of about one meter for best infrared pulse reception by the projector. While maintaining this position, hold the IN button (models CP 2033, CP 3000, and CP 3137) or the 1 button (models CP 2047, CP 2147, CP 1042, and CP 1142) down and with the other switch the projector on. Hold the remote button (IN or 1) down until the projector emits a sequence of short tones followed by one long tone signaling the end of the procedure. Release the remote button (IN or 1) and proceed with normal use.

Programming the remote control for use of multiple projectors in a single room.

Positioning the Projector and Focusing the Image on the Screen

For best quality when projecting the test images, and in particular the polarized tests, use the projection screen (17) supplied with the projector.

The screen must be placed at least 2.8 meters from the patient (not applicable in case of the R2033 Zoom projector) and used in dim lighting conditions or in any case in such a position as not to reflect sunlight (see also PROJECTION DISTANCE).

Once the projector and the screen are positioned, turn the ring (4) to focus the projector (models CP 2047, CP 2147, CP 1042, CP 1142, CP 3000, and CP 3137).

The procedure is different for the CP 2033 Zoom projector since this is the only model that permits conducting testing without having to align the patient with the projector (he may be up to 1 m in front of or behind the projector). For this reason, the user must be able to adjust the enlargement of the projected tests, besides the focus. To do so, open the top cover and adjust the two knobs: knob (25) to change the enlargement of the projected tests, and then (24) to focus the enlarged image. Once the desired settings are obtained, tighten the two knobs.



Image adjustment and focusing buttons – CP 2033 Zoom projector

WARNINGS, AND GENERAL SAFETY MEASURES

- Never touch the power cord with wet hands. Check frequently that the cord is so placed as not to be stepped on or crushed by weights. Never knot the cord.
- A damaged power cord can cause fires or electrical shocks. Check frequently that the instrument power cord is in good condition. If it becomes necessary to replace it, contact your supplier.
- Do not perform any repairs or maintenance work on the instrument or the electrical system beyond what is explained in this manual.
- Do not use the instrument near water and be careful not to spill liquids on any part of it. Avoid damp and dusty locations and locations subject to brusque changes in temperature and humidity.
- Disconnect the instrument from the line power supply before cleaning and/or disinfecting.
- **The instrument neither generates nor receives electromagnetic interference when used in proximity to other devices; no preventive or corrective measures need therefore be taken.**

List of remote programming buttons

Once the projector is correctly positioned and focused, proceed as follows to perform testing:

- a) Switch on the instrument with the main switch (7).
- b) Press a test button on the remote control to call up a test.
- c) Begin testing.

Attention! Never use the remote from a distance of less than 1.5 meters from the front of the remote. Lighting of a LED when a button is pressed on the remote indicates that the signal has been emitted and the batteries are charged. The projector will switch off automatically after 10 minutes of inactivity and will emit an audible signal 30 seconds prior to switch off (STANDBY).

Depending on projector model, the remote will contain different function buttons:

LAMP Button

When included, switches the projector lamp on or off. If the lamp is off, pressing any test button will cause it to switch on.

Programming Buttons

Each projector can store two programmed sequences of tests to be projected onto the screen, including masks (where present). The buttons used for programming are IN and OUT (CP 3000 and CP 3137), or the NEXT button and buttons 1 and 2. The programmed sequences, called up by pressing buttons 1 or 2, always begin with the first test stored in memory. To skip forward or backward press the NEXT and BACK buttons or the arrow buttons. The tests in a sequence may be substituted at any time.

Masking Buttons

All models except CP 1042 and CP 1142) permit using different masks for both horizontal and vertical lines and single letters or symbols.

Programming Test Sequences

To program test sequences on the remote, proceed as outlined below:

CP 3000 and CP 3137 Projectors

To memorize a sequence:

- 1) Press program button **1** or **2** (after programming, the same button will be used to call up the sequences).
- 2) Press **IN**.
- 3) Press the button for the first test to be memorized in the sequence.
- 4) Press the **NEXT** button once.
- 5) Repeat steps 3 and 4 until all the desired tests have been entered. After entering the last test, press the **NEXT** button once.
- 6) Press **OUT** to memorize the sequence on the button (**1** or **2**) selected in step 1.

To change a test in a programmed sequence:

- 1) Press program button **1** or **2**, according to the sequence to be modified.
- 2) Press the **NEXT** button to advance to the test to be substituted.
- 3) Press **IN**.
- 4) Press the button for the new test to be entered.
- 5) Press the **NEXT** button to advance to the end of the sequence. You will see an empty test circle.
- 6) Press **OUT**.

To interrupt a sequence:

- 1) Press program button **1** or **2**, according to the sequence to be modified.
- 2) Press **IN**.
- 3) Press the **NEXT** button to advance to the test at which the sequence is to be interrupted.
- 4) Press **OUT**.

CP 2047, CP 2147, CP 1042, CP 1142, and CP 2033 Projectors

To memorize a sequence:

- 1) With the projector on and in standby, hold down program button **1** or **2** for 5 seconds until a beep sounds to signal entry into programming mode. After programming, the same button will be used to call up the sequences.
- 2) Press the button for the first test to be memorized in the sequence.
- 3) Press **NEXT**.
- 4) Repeat steps 2 and 3 until all the desired tests have been entered.
- 5) After entering the last test, press **NEXT**.
- 6) Hold down the program button for 5 seconds until a beep sounds to signal that the projector has exited the programming mode and entered standby with the sequence memorized on button **1** or **2** selected in step 1.

To call up a sequence:

- 1) Press button **1** or **2**, according to the desired sequence.
- 2) The **NEXT** button advances to the next test in the sequence; the **BACK** button goes back to the previous test.
- 3) Press any key to exit the sequence.

To change a test in a previously-memorized sequence:

- 1) With the projector on and in standby, hold down program button for 5 seconds until a beep sounds to signal entry into programming mode.
- 2) Press **NEXT** to advance to the test to be substituted.
- 3) Press the button for the new test to be entered.
- 4) Press the **NEXT** button to advance to the end of the sequence. You will see an empty test circle.
- 5) Hold down program button **1** or **2** for 5 seconds until a beep sounds to signal that the projector has exited the programming mode and entered standby.


To interrupt a sequence:

- 1) With the projector on and in standby, hold down program button **1** or **2** for 5 seconds, according to the sequence to be modified.
- 2) Press the **NEXT** button to advance to the test at which the sequence is to be interrupted.
- 3) Hold down program button **1** or **2** for 5 seconds until a beep sounds to signal that the projector has exited the programming mode.


Projection Distance

The patient and the projector must be at the same distance from the projection surface, except in the case of the CP 2033 Zoom projector. The angle of vision, like measurement of visual acuity, is not influenced by this fact and the visus attributed each time will be valid for all distances. In order not to falsify the refraction value, the projection distance must fall between a minimum of 2.9 m and a maximum of 6.1 m.

With the Zoom projector (CP 2033) it is not necessary to align the projector and the patient with respect to the projection surface (tolerance is ± 1 m), but in this case the visus is dependent on the projection distance and verification on the projected enlargement will thus be required.

**Projection with Mirrors**

In rooms in which it is impossible to attain the minimum projection distance (2.9 m), the distance can be doubled by using mirrors.

The automatic chart projector with halogen lamp offers exceptionally clear images of even the smallest optotypes thanks to high-resolution optics and quality optotypes. The optotypes are sealed in a hermetic chamber to guarantee the cleanest possible test projections and stability over time.

The projector is controlled by an infrared remote unit in two modes, for selecting single tests with dedicated buttons or for programming sequences of tests.

The step/step motors are controlled by a microprocessor card for fast passage from one test to the next (0.03 seconds). The high-quality motors guarantee rapid, precise, vibration-free test presentation.

The halogen lamp switches off after 10 minutes of non-use (standby function).

CSO produces a wide range of projectors; below are presented some details of the different models, divided by category.

Variable-Contrast Chart Projector (CP 3000 and CP 3137)

These projectors permit varying contrast from 100% to 5% in linear progression with 10% steps to 10% and a final 5% step. This function permits administering special tests like comparison of visual acuity at high and low contrast and contrast sensitivity tests. The two-disk, 37-test CP 3000 projector features logMAR progression while the CP 3137 features decimal progression.

Dual-Disk Chart Projector, 30x enlargement (CP 2047 and CP 2147)

45-test projectors, red-green masking, polarized, horizontal-line, vertical-line, and single-letter tests. The CP 2047 projector features logMAR progression while the CP 2147 features decimal progression. The projectors can be interfaced, through the RS232 serial port, with the most common types of commercial computerized phoropters.

Single-Disk Chart Projector, 30x enlargement (CP 1042 and CP 1142)

42-test projector with no overlay filters or masks. The CP 1042 features logMAR progression while the CP 1142 features decimal progression. The projectors can be interfaced, through the RS232 serial port, with the most common types of commercial computerized phoropters.

Zoom Chart Projector (CP 2033)

33-test projector with red/green masks and polarized, horizontal- and vertical-line, and single letter tests. Progression according to Snellen notation (20 ft).

TECHNICAL SPECIFICATIONS

Technical Features	CP 2033	CP 2047, CP 2147, CP 1042, CP 1142	CP 3000, CP 3137
Projection distance	-	2,9/6,1 m	2,9/6,1 m
Tests	On optical glass disks		
Test and mask selection speed	0.03 seconds		
Projection enlargements	-	30 X at 5 m.	23,89 X at 5 m.
Tilt adjustment	Approx. ±. 15° in all directions		
Automatic switch off	10 minutes after last signal		
Projector dimensions	Width 205 mm – Length 340 mm – Height 210 mm		
Weight	5,5 Kg	5,5 Kg	6,5 Kg
Projector versions	For tabletop, wall, stand, or refraction unit installation		
Power supply	120Vac –230Vac 50-60 Hz		
Absorbed power	40 VA	40 VA	65 VA
Fuses	5x20 120V 500mAT 230V 250mAT		5x20 120V 1AT 230V 500mAT
Halogen lamps	12V 20W G4	12V 20W G4	2x 6V 25W G4

Classification According to EN60601-1 Standard

Class of protection against direct and indirect electrical shock:

Class I

Applied parts:

Type B

Compliance level for protection against humidity:

Common devices (shell not protected). IP20

Method of sterilization:

Disinfectable devices

Compliance level for protection in the presence of inflammable anesthetics or detergents:

No protection


Use conditions:

Continuous service


Patient-coupled electrical equipment definition:

Device with parts applied to patient


Block diagrams of CP 2047, CP 2147, CP 2033, CP 1042, and CP 1142 chart projectors.

Block diagrams of CP 3000 and CP 3137 chart projectors.

The manufacturer agrees to provide, on request, any diagrams, components lists, and/or detailed technical instructions that may be of use to authorized, previously-trained personnel for maintenance and calibration.


Lamp replacement on models CP 2047, CP 2147, CP1042, CP 1142, and CP 2033.

The projector requires no particular routine maintenance. The repairs illustrated below must all be performed only with the power cord disconnected from the line socket. For any type of breakdown requiring repairs different from those described below, contact your installation service.

Replacing the Lamp

Proceed as outlined below:

- Before beginning, disconnect the power cord from the line socket. Use the wrench (20) supplied as standard equipment to unscrew the two screws (5) holding the projector casing in place.
- Lift the casing (1).
- Lift and remove the lamp cover (11).
- Remove the unserviceable lamp. **Warning: the lamp may be very hot!**
- Replace with a new halogen lamp (10) or (12) in the case of the contrast projector, holding the lamp by the plastic envelope (NEVER touch the bulb with bare hands).
- Replace the lamp cover (11).
- Close the projector casing (1) and screw down (5).
- Reconnect the plug to the line socket.

Lamp replacement on models CP 3000 and CP 3137.

Replacing the Remote Control Batteries

Proceed as outlined below:

- Remove the cover on the back of the remote.
- Replace the batteries with 4 new 1.5V AAA batteries. The correct polarity is indicated on the inside of the battery compartment.
- Replace the cover.

Replacing the Line Fuses

To replace the line fuses:

- Before beginning, disconnect the power cord from the line socket.
- The line fuses are located at the rear of the instrument inside the socket/voltage selector unit (8).
- Extract the voltage selector and remove the unserviceable fuses.
- Replace the fuses with new fuses compatible with the line voltage supply as reported on the data plate on the back of the projector. Replace the voltage selector.
- Reconnect the plug to the line socket.


Protecting the Projector from Dust, Cleaning

When the instrument is not in use, it should be covered with the dust cover. Any dust that may accumulate on the projector should be periodically removed using a very soft cloth or a rubber ball air blower. Clean the outer surfaces with a cloth barely dampened with water. Do not use solvents.


Other Precautionary Measures

Unscheduled maintenance (repairs, replacement of components, checking internal parts, etc.) is the exclusive competence of the CSO technical assistance service.

Before sending the instrument to CSO for repairs or maintenance, clean all the surfaces and disinfect using a clean cloth dampened with ethyl alcohol.

Replacing the remote control batteries

Replacing the line fuses

CSO S.r.l. assumes responsibility for instrument compliance with European Directive 93/42/EEC and therefore for its performance, safety, and reliability.

CSO S.r.l. would like to point out that all the components of the devices it produces are covered by insurance and full warranty for 24 (twenty-four) months.

CSO S.r.l. nevertheless declines said responsibility if:

- installation and start-up are not performed in accordance with the instructions and precautions set forth in this manual.
- the instrument is used in ways not in accordance with the instructions and precautions set forth in this manual.
- accessories and/or spare parts not supplied or recommended by **CSO S.r.l.** are used..
- repairs and safety checks are not performed only by competent, qualified, and suitably-trained personnel authorized by **CSO S.r.l.**
- the electrical system of the instrument installation site does not comply with CEI standards and the pertinent laws and regulations in force.

CSO S.r.l. also declines any and all responsibility for direct or indirect consequences or damage to persons and/or things deriving from improper use of the instrument and/or from erroneous clinical evaluation of information derived from its use.

WARRANTY AND TECHNICAL ASSISTANCE

CSO S.r.l. guarantees this product for a period of 24 months from date of invoice. This warranty includes replacement, at CSO, of components and materials as well as relative labor. Shipping costs will be at the customer's expense.

Parts subject to wear and/or deterioration during normal use (for example, lamps and fuses) and parts damaged by improper use or inadequate maintenance are not covered by warranty.

CONDITIONS NOT INCLUDED UNDER WARRANTY

- Repair of damages caused by natural catastrophes, mechanical shock (dropping, crushing, etc.), defects in the user's electrical system, negligence, improper use, and/or maintenance/repairs performed using non-original materials and/or by persons not authorized by **CSO S.r.l.**
- Any type of improper use or use not specifically intended by the manufacturer.


CSO S.r.l. declines responsibility for any interruption or inefficiency in service due to causes or circumstances beyond its control. The customer shall in no case have any right to compensation for damages suffered as a consequence of the unavailability of the instrument.

To request technical assistance with maintenance, contact a technical assistance center or directly contact:

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Fuse

Applied parts classified as Type B in accordance with EN 60601-1 standard

General warning of the duty to read the instruction manual carefully before installing and using the instrument.

"CE marking" attesting to product compliance with European Union Directive 93/42/EEC ("Medical Devices").

C.S.O. srl - Costruzione Strumenti Oftalmici with headquarters in V. Degli Stagnacci 12/e - 50010 Badia a Settimo - Firenze – Italia, manufacturers of the medical electrical device (Class I, as per Directive 93/42/EEC) denominated

CHART PROJECTOR

MODEL: CP 2047, CP 2147, CP 1042, CP 1142, CP 2033, CP 3000, CP 3137,

in the person of its legal representatives Sergio Mura and Giuseppe Matteuzzi, who assume full personal liability for the following,

WHO ASSUMES

full personal liability for the following, hereby certifies

- complies fully with the following EU directive, including the latest amendments, and with the relative national implementing legislation:

DIRECTIVE 93/42/EEC "MEDICAL DEVICES " of 14 June 1993

and that

- the following harmonized standards were applied:

- EN 60601-1 - "MEDICAL ELECTRICAL EQUIPMENT – PART 1: GENERAL REQUIREMENTS FOR SAFETY" as amended, corresponding to Italian standard CEI 62-5 - 2005 edition
- EN 60601-1-2 - "MEDICAL ELECTRICAL EQUIPMENT – PART 2: Collateral Standard: Electromagnetic Compatibility," corresponding to Italian standard CEI 62-50 – 2001 edition;
- UNI EN ISO 8596 (1996) "Visual Acuity Test Types. Standard Optotype and its Presentation."
- UNI EN ISO 8597 (1996) "Visual Acuity Test Types. Method for Correlating Optotypes."
- ISO 10938: (1998) "Ophthalmic Instruments. Chart Projectors."

- The instrument is marketed with the  mark attesting to compliance with the above-mentioned directives.

SCANDICCI (FIRENZE), 28/09/06

The legal representatives of CSO Srl
(Sergio Mura and Giuseppe Matteuzzi)

.....

The following reference standards were applied in design, production, and testing of the product:

EU Directives

- **DIRECTIVE 93/42/EEC "MEDICAL DEVICES " of 14 June 1993.**
- **DIRECTIVE 2002/96/EC "Waste Electrical and Electronic Equipment (WEEE)."**

Standards concerning Quality Management Systems

- **UNI EN ISO 9001:2000 "Quality Management Systems – Requirements."**
- **UNI EN ISO 13485:2004 "Medical Devices – Quality Management Systems - Clinical Requirements for Regulatory Compliance."**

Technical Standards

- **EN 60601-1 - "MEDICAL ELECTRICAL EQUIPMENT – PART 1: GENERAL REQUIREMENTS FOR SAFETY" as amended, corresponding to Italian standard CEI 62-5 - 2005 edition.**
- **EN 60601-1-2 - "MEDICAL ELECTRICAL EQUIPMENT – PART 2: Collateral Standard: Electromagnetic Compatibility," corresponding to Italian standard CEI 62-50 – 2001 edition.**
- **UNI EN ISO 8596 (1996) "Visual Acuity Test Types. Standard Optotype and its Presentation."**
- **UNI EN ISO 8597 (1996) "Visual Acuity Test Types. Method for Correlating Optotypes."**
- **ISO 10938: (1998) "Ophthalmic instruments. Chart Projectors"**
- **UNI EN ISO 15004 "Ophthalmic instruments. Fundamental Requirements and Test Methods." 2000 edition.**
- **UNI EN ISO 14971:2004 "Risk Management for Medical Devices."**

The registered trademarks are property of the respective trademark owners.

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CSO Srl is certified according to